


The Inflatable  
Lumbo-Sacral  
Orthoses (I.L.S.O.)  
from PMT.

INFLATABLE LUMBO-SACRAL ORTHOSES

**I.L.S.O.**



PMT specializes in  
designing custom  
ILSO's to fit your  
patient's needs.

### *Benefits...*

- Moveable posterior airbag ensures customized spinal support.
- Easy access valve gives patient control of desired airbag pressure.
- Lightweight, low-profile design.
- Bi-valve construction for greater A-P compression.
- Quick and easy fitting.
- Pendulous anterior shells available.

**PMT**<sup>®</sup>

*Quality, Innovation, Customer Driven*

# Product Information

# I.L.S.O.

## Description:

### Model 1515 PMT ILSO/Inflatable Lumbo-Sacral Orthosis (I.L.S.O.)

Complete System includes - Anterior and posterior shells, airbag, valve manifold, straps and inflation bulb.

### Model 1525 PMT IITLSO/Inflatable Thoracic Lumbo-Sacral Orthosis

Complete System includes - Anterior and posterior shells, sternal pad, shoulder straps, airbag, valve manifold, straps and inflation bulb.

### Model 1530 PMT IITLSO/Inflatable Thoracic Lumbo-Sacral Orthosis

Complete System includes - Anterior and posterior shells, airbag, valve manifold, straps and inflation bulb.

## Sizes:

PMT stocks standard male and female sizes. Including pendulous anterior shells. PMT's sizing chart is available upon request.

## Measurements Needed:

- Chest circumference at xiphoid process
- Waist circumference
- Hip circumference
- Xiphoid to pubis (vertical measurement)
- Sternal notch to pubis (for Model 1525 & 1530)

## Note:

Custom braces available upon request. All customs shipped within 2 business days of receiving measurements.

▶ **Local sales representatives are available to serve you. Contact PMT Corporation for more information.**

PMT and MANKIND are registered trademarks of PMT Corporation  
Copyright © 2004. All rights reserved.

**PMT**<sup>®</sup>  
*Quality. Innovation. Customer Driven*

## ILSO/ITLSO Model Number Matrix

System	Without Foam Lining	With Foam Lining
2 Strap System	1515-1 ILSO 1525-1 IITLSO w sternal pad 1530-1 IITLSO	1515-2 ILSO 1525-2 IITLSO w sternal pad 1530-2 IITLSO
3 Strap System	1515-3 ILSO 1525-3 IITLSO w sternal pad 1530-3 IITLSO	1515-4 ILSO 1525-4 IITLSO w sternal pad 1530-4 IITLSO

For assistance with measurements, contact your PMT representative. PMT is willing to customize any ILSO/ITLSO to fit your needs.



Model = 1515-4 ILSO on left and Model =1530-4 IITLSO shown at right.

## Our Commitment... Customer Satisfaction

PMT is committed to providing you the highest quality products and the finest customer service. We are available to assist you Toll-Free 24 hours a day, 7 days a week, 365 days a year. If it's more convenient, you can even fax us your requirements.

PMT CORPORATION  
1500 Park Road  
Chanhassen, MN USA 55317

Call Toll Free... 24 Hours a Day/365 Days a Year  
**1-800-MANKIND**<sup>®</sup> (626-5463)  
Phone 952-470-0866 Fax 952-470-0865