


INFLATABLE THORACIC LUMBO-SACRAL ORTHOSIS

I.T.L.S.O. MODEL 1535

- New innovative gear design
- Adjustable or anchored strap construction available
- Easy access valve gives patient control of desired airbag pressure
- Light weight low profile design
- Bi valve construction for greater A-P compression
- Pendulous anterior shells available
- Adjustable thoracic assembly and high back posterior shell provide greater patient comfort and support



- Innovative gear design allows for ease of adjustment



- Adjustable thoracic bar provides greater comfort and support



- Moveable posterior airbag ensures customized spinal support



- Quick and easy fitting

PMT®

Quality, Innovation, Customer Driven

Product Information

MODEL 1535
I.T.L.S.O.

Description:

Model 1535 PMT ITLSO/Inflatable Thoracic
Lumbo-Sacral Orthosis (I.T.L.S.O.)

Complete system includes - anterior and posterior shells, thoracic bar with adjustable gear attachment, air bag valve manifold, straps and inflation bulb.

1535-1 Available with adjustable straps

1535-2 Available with anchored straps

Sizes:

PMT stocks standard male and female sizes. Including pendulous anterior shells. PMT's sizing chart is available upon request.

Measurements Needed:

- Chest circumference at xyphoid process
- Waist circumference
- Hip circumference
- Xyphoid to pubis (vertical measurement)
- Sternal notch to pubis

Note:

Custom braces available upon request.

All customs shipped within 2 business days of receiving measurements.

▶ Local sales representatives are available to serve you. Contact PMT Corporation for more information.

Our Commitment... Customer Satisfaction

PMT is committed to providing you the highest quality products and the finest customer service. We are available to assist you Toll-Free 24 hours a day, 7 days a week, 365 days a year. If it's more convenient, you can even fax us your requirements.

PMT and MANKIND are registered trademarks of PMT Corporation.
© Copyright 2007 PMT Corporation. All Rights Reserved

PMT®

Quality, Innovation, Customer Driven

PMT CORPORATION
1500 Park Road
Chanhassen, MN USA 55317
www.pmtcorp.com

Call Toll Free... 24 Hours a Day/365 Day a Year
1-800-MANKIND® (626-5463)
Phone 952-470-0866 Fax 952-470-0865