

OSADA Dermabrasion System (Model XL-30W)

- **Power Control Box with handle** - lightweight and portable
 - simple to use controls
 - on/off and forward/reverse switches
 - manual setting for constant speed
 - wide speed range (1,500 to 30,000 RPM)
 - international voltage switch (120V or 220V)
- **Autoclavable Handpiece** - smooth and vibration free
 - protective nose guard
 - lever chuck for easy bur changing and locking
 - accepts all standard size bits, fraises and scalp transplant punches
- **Variable speed foot pedal**
- **Six foot coiled motor cord**



Distributed by **PMT**® / **Permark**®

1500 Park Road, Chanhassen, Minnesota

952-470-0866 Telephone • 952-470-0865 Facsimile • 800-PERMARK Toll Free • www.permark.com

©2005 PMT Corporation. All Rights Reserved.