


## Needles

Sterile. Sharp, smooth, tapered points. 12 different needles for cosmetic and reconstructive procedures. Made from corrosion-resistant stainless steel. Lot coded. Autoclavable.

## Needle Suggestions

Size	Description	Suggested Uses
#1	Single need	Fine eyebrow, lip liner, small area camouflage
#3 Tight	Three fine needles in a cluster	Fine eye liner, lip liner, small area camouflage
#3	Three needles in a cluster	Eye liner, lip liner, small area camouflage
#4 Flat	Four needles in a single row	Eyebrow, mall area camouflage, shading
#5	Five needles in a cluster	Heavy eye liner, lip liner, camouflage
#6 Flat	Six needles in a single row	Nipple, areola, camouflage, shading
#7	Seven needles in a cluster	Nipple, areola, heavy eye liner, lip liner, camouflage
#7 Flat	Three needles in a single row + Four needles in a single row	Areola, camouflage, shading
#9	Nine needles in a cluster	Nipple, areola, camouflage
#11 Flat	Five needles in a single row + Six needles in a single row	Areola, camouflage, shading
#14	Fourteen needles in a cluster	Areola, camouflage
#18	Eighteen needles in a cluster	Areola, camouflage



**Permark**<sup>®</sup>  
Division of **FMT**

**COLOR CHART**  
FOR MICROPIGMENTATION

Contact us at 1-800-MANKIND, [www.permark.com](http://www.permark.com) or by fax 952-470-0865.  
PMT Corporation. 1500 Park Road, Chanhassen, MN 55317 USA



## Color Suggestions

### Pigments

Sterile, thick, creamy concentrated formula. 52 shades for cosmetic and reconstructive procedures. No water. No splattering. No mess. Made with purified, FD&C ingredients. Hypoallergenic. Lot coded. Dated for freshness. Quality made in the USA since 1985.

### Warning

While adverse reaction to pigment is rare, patch testing is recommended. Colors are not always permanent. They may change or fade when placed in skin. Black, Dark Gray, Navy Blue and Slate Blue are not recommended for eyebrow procedures. Terra Cotta, Plum and Blush have been reported to appear blue in the lips.

### Camouflage

- Flesh .5
- Flesh 1
- Flesh 1.5
- Flesh 2
- Flesh 2.5
- Flesh 4
- Flesh 9
- Flesh 10
- Medium Brown
- Sandstone 1
- Sandstone 2
- Sandstone 3

### Nipple-areola

- Dark Brown
- Flesh 2.5
- Flesh 4
- Flesh 5
- Flesh 6
- Flesh 7
- Flesh 8
- Flesh 9
- Flesh 10
- Flesh 11
- Flesh 12
- Light Brown
- Medium Brown
- Topaz
- Umber

### Eye Liner

- Black
- Brown Black
- Charcoal Green
- Dark Gray
- Green
- Medium Gray
- Navy Blue
- Pewter
- State Blue
- Umber

### Eyebrow

- Blonde
- Blonde K21
- Chocolate
- Dark Chocolate
- Kajal
- Taupe

### Lip Colors

- Autumn Red
- Blush
- Brick
- Coral
- Desert Rose
- Mauve Mocha
- Night Red
- Plum
- Red
- Terra Cotta

### Correctional

- Blonde
- Brown Toner
- Flesh .5
- Red Toner
- White Toner
- Yellow Toner