


NEW PMT[®] SURGICAL SUCTION INSTRUMENTS



CONICAL DESIGN WITH INTERNAL TAPERING

- » Allows greater visibility of the surgical site
- » Allows evacuation of debris through the entire suction site
- » Reduces obstructions from clogging the instrument
- » Includes a coiled style that matches taper

A TRULY MALLEABLE INSTRUMENT

- » High grade German stainless steel
- » Extra long life span

60 DIFFERENT SIZES

- » 8 different French sizes 3 to 10
- » 6 different lengths from 80mm to 175mm

EASY TO UNDERSTAND NUMBERING SYSTEM

- » French sizes and length of the instrument are clearly printed on the suction control chamber

TEARDROP SHAPED SUCTION CONTROL

- » For added comfort of the physician during long hours of use
- » Precise finger control of suction for maximum sensitivity

PRODUCT INFORMATION

NEW PMT[®] SURGICAL SUCTION INSTRUMENTS



Order Information:

To order PMT[®] Surgical Suction Instruments use the following table to determine the model number needed.

9001 - Micro suckers with French sizes 3–6

9002 - Macro suckers with French sizes 7–10

Model# Example 9001-04-100 (micro sucker French 4, 100mm length)

French

Length	3	4	5	6	7	8	9	10
80mm	9001-03-80	9001-04-80	9001-05-80	9001-06-80	9002-07-80	9002-08-80	9002-08-80	9002-10-80
100mm	9001-03-100	9001-04-100	9001-05-100	9001-06-100	9002-07-100	9002-08-100	9002-09-100	9002-10-100
120mm	9001-03-120	9001-04-120	9001-05-120	9001-06-120	9002-07-120	9002-08-120	9002-09-120	9002-10-120
145mm	9001-03-145	9001-04-145	9001-05-145	9001-06-145	9002-07-145	9002-08-145	9002-09-145	9002-10-145
160mm	9001-03-160	9001-04-160	9001-05-160	9001-06-160	9002-07-160	9002-08-160	9002-09-160	9002-10-160
175mm	9001-03-175	9001-04-175	9001-05-175	9001-06-175	9002-07-175	9002-08-175	9002-09-175	9002-10-175



Local sales representatives are available 24 hours a day, 7 days a week, 365 days a year to assist you. Contact PMT[®] Corporation for more information.

FDA REGISTERED

© Copyright 2014 PMT[®] Corporation. All Rights Reserved

PMT[®]
PMT CORPORATION
QUALITY, INNOVATION AND CUSTOMER DRIVEN

PMT CORPORATION
1500 Park Road
Chanhassen, MN USA 55317
www.pmtcorp.com

Call Toll Free ... 24 Hours a Day/365 Day a Year
1-800-MANKIND[®] (626-5463)
Phone 952-470-0866 » Fax 952-470-0865

PN. 22101 Rev. 1, 03/14